



Debrief from the roadside to the bedside

SimView™ Mobile

Increase the educational impact of debriefing your in-situ simulations by incorporating synchronized audio, video, and simulation information. SimView™ Mobile is the portable resource for an integrated capture and debriefing system. SimView Mobile can be easily transported and set up anywhere you run simulations to record the audio, 3 video streams, and a patient monitor stream in one debriefing file.

SimView™ Mobile



Capture

- Grab-n-Go with SimView Mobile's backpack to record in varied environments
- Designed for portability and simplicity
- Quickly and easily set up the cameras and audio
- Run LLEAP software from SimView Mobile to control simulations while you record



Share

- Share files with faculty and students
- Free Session Viewer software allows staff and students to playback files and edit/annotate event log
- Upload files from SimView Mobile to SimView server for storage and or playback
- Integrate debrief files with SimManager



Debrief

- Debrief onsite immediately after your in-situ simulation or save the file to debrief later
- Annotate event logs for comments, medications, events, correct/incorrect steps
- Move back and forth along the timeline to highlight specific areas for review



Implementation and support

SimView Mobile is pre-configured to run when you are ready to implement. From the included Quick Set-Up Guide, you simply match color coded connections established for the cameras and hub to form the system. Your audio source is also pre-set so you can position it where you want to pick up interactions.

You can also utilize SimView Mobile as your instructor application. Run your Laerdal simulator with the included LLEAP software and also pull in up to 3 video streams and a patient monitor.

Complement your SimView server

SimView Mobile is an excellent companion product to SimView server. SimView Mobile allows you the flexibility to move out of your simulation lab to run in-situ simulations, whereas SimView server is traditionally used for fixed simulation centers. Your recordings from SimView Mobile can be uploaded to SimView server.

Ordering Information:

216-20001 SimView Mobile

SimView Mobile includes: SimView Mobile software, laptop, 3 cameras, camera grips, USB HUB, external speaker, cabling, and backpack.

Optional Accessories include:

216-06010 SimView Mobile portable projector

216-06040 SimView Mobile Tripod

216-06050 SimView Mobile Hard case